

Sleep masks

Our complete line of nasal, full face, minimal contact and pediatric masks reflects our commitment to innovation and meeting patient and business needs.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



DreamWear

Like wearing nothing at all¹⁻³

Patients can sleep how they want with the freedom and comfort of DreamWear

The Philips Respironics DreamWear frame supports full face, nasal and gel pillows cushions to meet patients' needs. DreamWear's remarkable design directs airflow through the frame so patients can sleep comfortably, and the tubing connection located at the top of the head allows patients to sleep in any position they desire – on their stomach, side or back.



Users told us that DreamWear gave them more freedom to move, had a better fit than their prescribed mask and comes closest to making them feel like they don't have anything on their face during therapy.¹⁻³

Give your patients the headgear with all-night staying power

Our innovative, slip-resistant design – similar to eyeglasses holding the mask securely in place above, but not touching the ears – is designed to provide more stability and an improved fit and seal without compromising the comfort that your patients expect from the DreamWear mask.



In a recent study, 98% of gel pillows⁴ and 94% of nasal users⁵ stated that the updated headgear on DreamWear masks stays in place during use.

Note: Switching from a nasal cushion to a full face cushion requires different headgear and instructions. Consumers must consult their provider before making adjustments.
¹DreamWear Under the Nose Nasal. Data analysis at 30 days of use of 2015 Philips sponsored patient preference trial (n=98). Prescribed masks include ResMed Swift FX, ResMed Mirage FX, Philips Wisp, and ResMed P10. Data analysis at 30 days of use of 2019 Philips sponsored patient preference trial (n=111). Prescribed masks include ResMed AirFit N10, N20 Mirage FX, Philips Wisp, Pico, ComfortGel Blue, F&P Esor 5, Esor 2.
²DreamWear Full Face. Data analysis of Aug 2017 Patient Preference trial where (n=85) and prescribed masks include ResMed AirFit P10 and P20, Respironics Airova Now and F&P Simplicity.
³DreamWear Gel Pillows. Data analysis at 30 days of use of 2017 Philips sponsored patient preference trial (n=80). Prescribed masks include ResMed Swift FX and P10, F&P Opus and P420, and Philips Nuance and Nuance Pro. Data analysis at 30 days of use of 2019 Philips sponsored patient preference trial (n=72). Prescribed masks include ResMed AirFit P10, Swift FX, Philips Nuance / Nuance Pro, F&P Brovada, Opus 5, P10 Pro.
⁴Robot Preference Study with DreamWear gel pillows with original headgear users. (n=47), Jan/March 2019, at 30 days. Original headgear references to the sit headgear design on DreamWear gel pillows mask.
⁵Robot Preference study with DreamWear Nasal with original headgear users. (n=40), Jan/March 2019, at 30 days. Original headgear references to the sit headgear design on DreamWear UTN Nasal mask.

One mask

fits all three cushion types

Simply switch cushions without changing masks

Choosing the right mask for your patients is one of the most difficult jobs you face. The non-stop trial and error involving the masks, cushions and headgear needed to fit your patients becomes tedious and time-consuming.

On top of that, you must manage patients in their initial acceptance period who struggle for a variety of reasons including poor fit, leaking, nasal congestion and others.

Now, DreamWear's design allows you to quickly and easily switch between full face, nasal and gel pillows cushions without changing masks. With fewer mask decisions to make, you can focus on what's really important: your patients.

Benefits

- Can accommodate full face, nasal and gel pillows cushions
- Allows for easy switching between cushion types
- May reduce the number of masks needed to set up a patient
- Helps streamline the amount of inventory to be managed
- Can increase efficiency by requiring education on only one mask




Clinicians rated the DreamWear mask with nasal and gel pillows cushions as taking less time to change than another pillows mask and nasal mask.*

*Data analysis of RF and ST (n=20) survey responses regarding the amount of time needed to switch between two different masks on a patient.



DreamWear



Soft silicone frame
Light, flexible material designed for comfort.

Tip-of-head tube connection
Provides freedom of movement throughout the night.

Open field of vision
Wear glasses, read and watch TV.

In-frame airflow
Unique design allows airflow through the frame to help patients sleep comfortably.

No red marks*
Leaves no marks on the bridge of the nose.

Headgear with arms
Slip-resistant design provides more stability.

Innovative
Under-the-nose design.



Interchangeable cushions
Easily switch between full face, nasal and gel pillows cushions.

*Mask does not directly contact the bridge of the nose or nostrils.



	Nasal part number	Gel pillows part number	Full face part number
DreamWear ftpack			
DreamWear with headgear, medium frame, all cushion sizes	1116700	1124984	1133400
DreamWear without headgear, medium frame, all cushion sizes	1116701	1124985	N/A
DreamWear mask without headgear, medium frame			
Small cushion	1116705	1124986	1133405
Medium cushion	1116706	1124987	1133406
Large cushion	1116707	1125008	1133407
Medium-wide cushion*	1116708	N/A	1133408
DreamWear mask without headgear, small frame			
Small cushion	1116710	1125009	1133410
Medium cushion	1116711	1125010	1133411
Large cushion	1116712	1125011	1133412
Medium-wide cushion*	1116713	N/A	1133413
DreamWear mask without headgear, large frame			
Small cushion	1116715	1125012	
Medium cushion	1116716	1125013	
Large cushion	1116717	1125014	
Medium-wide cushion*	1116718	N/A	
DreamWear mask with headgear, small frame			
Small cushion	1116685	1125018	1133385
Medium cushion	1116686	1125019	1133386
Large cushion	1116687	1125020	1133387
Medium-wide cushion*	1116688	N/A	1133388
DreamWear mask with headgear, medium frame			
Small cushion	1116680	1125015	1133380
Medium cushion	1116681	1125016	1133381
Large cushion	1116682	1125017	1133382
Medium-wide cushion*	1116683	N/A	1133383
DreamWear mask with headgear, large frame			
Small cushion	1116690	1125021	1133390
Medium cushion	1116691	1125022	1133391
Large cushion	1116692	1125023	1133392
Medium-wide cushion*	1116693	N/A	1133393
Replacement parts			
Small cushion	1116740	1125035	1133430
Medium cushion	1116741	1125036	1133431
Large cushion	1116742	1125037	1133432
Medium-wide cushion*	1116743	N/A	1133433
Small frame	1116745	1116745	1116745
Medium frame	1116746	1116746	1116746
Large frame	1116747	1116747	1116747
Headgear	1116750	1116750	1133450
Fabric wraps	1116754	1116754	1116754
Sizing gauge	1116752	1125038	1116752

*A medium-wide cushion is not available for gel pillows.



DreamWisp

Comfort with every turn

Introducing DreamWisp minimal contact mask from Philips Respironics, the #1 sleep therapy* brand¹ of U.S. sleep physicians. DreamWisp helps your patients sleep comfortably with every turn because it marries the freedom of DreamWear's revolutionary, top-of-the-head tube design with the seal of Wisp's minimal contact nasal cushion.

Now, you can provide your patients with an innovative solution that lets them sleep how they want without having restrictive tubing on the face. DreamWisp's design directs airflow through the frame so patients can sleep comfortably, and the tubing connection at the top of the head allows patients to sleep in any position they want – on their stomach, side or back.

*CPAP and masks: Philips Respironics

BB 800.345.6443 www.philips.com/respironics

[Go to table of contents](#)

Benefits of DreamWisp

- Tubing at the top-of-the-head allows patients to sleep in any position they want.
- Airflow through the mask frame allows patients to sleep comfortably
- Minimal contact nasal cushion provides comfort and seal
- Four-point headgear provides stability

DreamWisp nasal mask with medium connector fitpak configurations

	Part number
DreamWisp with headgear, small, medium and large cushions	1137916
DreamWisp without headgear, small, medium and large cushions	1137918

DreamWisp nasal mask without headgear with medium connector

Petite DreamWisp cushion	1137922
Small DreamWisp cushion	1137923
Medium DreamWisp cushion	1137924
Large DreamWisp cushion	1137925
Extra-large DreamWisp cushion	1137926

DreamWisp nasal mask with headgear with medium connector

Petite DreamWisp cushion	1137932
Small DreamWisp cushion	1137933
Medium DreamWisp cushion	1137934
Large DreamWisp cushion	1137935
Extra-large DreamWisp cushion	1137936

DreamWisp nasal mask cushions

Petite DreamWisp cushion	1137956
Small DreamWisp cushion	1137957
Medium DreamWisp cushion	1137958
Large DreamWisp cushion	1137959
Extra-large DreamWisp cushion	1137960

Replacement parts

DreamWisp small connector	1137961
DreamWisp medium connector	1137962
DreamWisp large connector	1137963
DreamWisp headgear	1137964
DreamWisp large headgear	1137965
DreamWisp RS headgear	1137966
DreamWisp sizing gauge	1137967
DreamWisp fabric wraps	1137968
DreamWisp magnetic clips	1120434
DreamWisp frame faceplate	1137970
DreamWisp EE elbow	1116748



Available in 5 sizes:



Order online, anytime at www.homehealthstore.philips.com

800.345.6443 89

[Go to table of contents](#)



Wisp

A combination of comfort, performance and style

Developing solutions that meet the needs of a rapidly changing industry calls for fresh thinking and collaboration. At Philips Respironics, we traveled around the world asking clinicians, homecare providers and patients how we could improve their experience.

Clinicians wanted us to resolve the tradeoff between nasal and pillows masks, homecare providers wanted a solution that could reduce operational complexity by getting more patients served from one package and patients wanted a better overall mask experience that improved their lifestyle.

Key benefits

- Compact design for a comfortable fit
- Soft and comfortable for a restful night's sleep
- Open field of vision for a feeling of independence
- Quick and easy to adjust for a superior seal and personalized fit
- Contemporary frame options to complement your style

90 800.345.6443 www.philips.com/respironics

[Go to table of contents](#)

Part number	HCPCS	
Wisp, fabric frame		
1094051	Wisp mask with fabric frame, with headgear	A7034, A7035
1094056	Wisp mask with fabric frame, without headgear	A7034
Wisp, clear frame		
1094050	Wisp mask with clear frame, with headgear	A7034, A7035
1094055	Wisp mask with clear frame, without headgear	A7034
Wisp youth		
1109298	Wisp youth fabric, P, S/M, L cushions with reduced size headgear	
Wisp pediatric		
1104953	Fit pack with 3 cushion sizes	
Wisp pediatric replacement parts		
1104973	Headgear	
1104977	Elbow/tube with cover	
1104972	Frame	
1104969	Small cushions (SCS)	
1104970	Medium cushions (SCM)	
1104971	Large cushions (SCL)	
Wisp replacement parts		
1094086	Petite cushion	A7033
1094087	Small/medium cushion	A7033
1094088	Large cushion	A7033
112031	Extra-large cushion	
1094080	Clear frame	
1094081	Fabric frame	
1094082	Wisp headgear	A7035
1094078	Wisp headgear, large size	A7035
1109307	Wisp headgear, reduced size	
1094083	Wisp headgear clips	
11200434	Wisp magnetic headgear clips	



Adult: All three sizes in one package



Youth: All three sizes in one package



Compact design comfortably fits more than **98%** of patients



Open field of vision for a feeling of **greater independence** – patients can read and even wear their glasses

Order online, any time at www.homehealthestore.philips.com

800.345.6443 91

[Go to table of contents](#)



Pico

Pico nasal mask delivers performance for you and your patients

Practical innovation

At Philips Respironics, we created Pico with a definitive goal in mind – to create a truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal. These qualities make it an ideal choice for new and existing sleep apnea patients who want minimal disruptions to their lifestyle.



Part number		HCPCS
Pico		
1104940	Pico nasal mask fitpack, small/medium, large, extra-large, with headgear	A7034, A7035
1104915	Pico nasal mask, small/medium, with headgear	A7034, A7035
1104916	Pico nasal mask, large, with headgear	A7034, A7035
1104917	Pico nasal mask, extra-large, with headgear	A7034, A7035
1104918	Pico nasal mask, small/medium, without headgear	A7034
1104919	Pico nasal mask, large, without headgear	A7034
1104920	Pico nasal mask, extra-large, without headgear	A7034
Replacement parts		
1104934	Pico headgear	A7035
1104935	Pico headgear, reduced size	A7035
1104936	Small/medium nasal cushion	A7032
1104937	Large nasal cushion	A7032
1104938	Extra-large nasal cushion	A7032
1104939	Pico sizing gauge	A7032
1104933	Pico headgear clips	
1104931	Pico quick release swivel	

Available in fitpack and sizes: **S/M** **L** **XL**
 92 800.345.6443 www.philips.com/respironics

Order online, anytime at www.homehealthstore.philips.com 800.345.6443 93

[Go to table of contents](#)



ComfortGel Blue



ComfortGel Blue

Staying true to our acclaimed blue gel masks, we've retained the globally popular comfort and flexibility features you expect, while enhancing patient comfort even further. With ComfortGel Blue, we've reduced noise, and redirected exhalation air flow away from the bed partner. These and other enhancements may help to create a better patient experience and encourage greater adherence.

Part number		HCP5
ComfortGel Blue with headgear		
1070040	Petite	A7034, A7035
1070039	Small	A7034, A7035
1070038	Medium	A7034, A7035
1070037	Large	A7034, A7035
ComfortGel Blue without headgear		
1070044	Petite	A7034
1070043	Small	A7034
1070042	Medium	A7034
1070041	Large	A7034
ComfortGel Blue with headgear DuoPack (includes two same-size cushion sets)		
1070053	Small and medium	A7034, A7035, A7032
ComfortGel Blue replacement parts		
1070108	Petite cushion and flap	A7032
1070107	Small cushion and flap	A7032
1070106	Medium cushion and flap	A7032
1070105	Large cushion and flap	A7032
1070104	Petite cushion	A7032
1070103	Small cushion	A7032
1070102	Medium cushion	A7032
1070101	Large cushion	A7032
1033678	Premium headgear	A7035
1059236	RS Premium headgear (smaller)	A7035
1040114	Forehead pad	
1009723	ComfortGel Blue Quick Clips (2 per pack)	

Amara View

Enjoy the view

Amara View has a wider field of vision than Resmed AirFit F20 and Fisher & Paykel Simplus.**

An idea whose time has come. Amara View's innovative design prevents red marks, discomfort or irritation on the nose bridge.* Amara View is smaller, lighter†, and covers less of a patient's face than Resmed AirFit F20 and Fisher & Paykel Simplus.** It's easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Key benefits

- Innovative design prevents red marks, discomfort or irritation on the nose bridge
- Enables wearing glasses, reading, and watching TV
- Easy to use with one modular frame for all cushion sizes
- No bulky cushion and frame in front of your face and eyes



**Design claim
***In-lab testing (October 2013, November 2013, May 2014), comparing Amara View with Quattro, Quattro FX, Quattro Air, APPE F30 and F&P Simplus. †Without open release tube.

Part number

HCPCS

Part number	Description	HCPCS
1090622	Amara View small mask with headgear	A7030, A7035
1090623	Amara View medium mask with headgear	A7030, A7035
1090624	Amara View large mask with headgear	A7030, A7035
1090632	Amara View small mask without headgear	A7030
1090633	Amara View medium mask without headgear	A7030
1090634	Amara View large mask without headgear	A7030
1090670	Amara View Fit Pack small, medium and large cushions with headgear	A7030, A7035

Amara View replacement parts

1090689	Amara View quick-release tube	
1090692	Amara View small cushion	A7031
1090693	Amara View medium cushion	A7031
1090694	Amara View large cushion	A7031
1090696	Amara View headgear, large	A7035
1090697	Amara View headgear, standard	A7035
1090698	Amara View headgear clip, 2 pk	
1090699	Amara View sizing gauge	
1136626	Amara View magnetic clip, 2 pk	

96 800.345.6443 www.philips.com/respirators

[Go to table of contents](#)



Order online, anytime at www.homehealthstore.philips.com

800.345.6443 97



Amara

Performance, appeal and simplicity at its best

With gel and silicone options and our one-“click” design, Amara delivers performance and ease of use for your patients. With the addition of a gel cushion option, Philips Respironics has made the ingeniously-simple Amara full face mask even more brilliant. Your patients can now choose the type of comfort they prefer in their cushion with gel or silicone options.

Silicone

Part number		HCPCS
Amara with headgear		
1090200	Petite (reduced size frame and headgear)	A7030, A7035
1090201	Small (reduced size frame and headgear)	A7030, A7035
1090226	Medium (reduced size frame)	A7030, A7035
1090227	Large (reduced size frame)	A7030, A7035
1090267	Amara silicone fit pack (includes all four silicone cushions, reduced size frame and headgear)	
Amara without headgear		
1090210	Petite (reduced size frame)	A7030
1090211	Small (reduced size frame)	A7030
1090215	Medium (reduced size frame)	A7030
1090216	Large (reduced size frame)	A7030
Amara replacement parts		
1090290	Petite cushion	A7031
1090292	Small cushion	A7031
1090293	Medium cushion	A7031
1090294	Large cushion	A7031
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1078777	Talon headgear clips	
1081899	Entrainment valve swivel	



Available in 4 sizes:



Gel

Part number		HCPCS
Amara gel with headgear		
1090400	Petite (reduced size frame and headgear)	A7030, A7035
1090401	Small (reduced size frame and headgear)	A7030, A7035
1090405	Medium (reduced size frame)	A7030, A7035
1090406	Large (reduced size frame)	A7030, A7035
1090467	Amara gel fit pack (includes all four gel cushions, reduced size frame and headgear)	
Amara gel without headgear		
1090410	Petite (reduced size frame)	A7030
1090411	Small (reduced size frame)	A7030
1090415	Medium (reduced size frame)	A7030
1090416	Large (reduced size frame)	A7030
Amara gel replacement parts		
1090490	Petite cushion	A7031
1090492	Small cushion	A7031
1090493	Medium cushion	A7031
1090494	Large cushion	A7031
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1078777	Talon headgear clips	
1081899	Entrainment valve swivel	
Amara starter kits (includes gel mask with same size silicone cushion)		
1106618	Petite	A7030, A7035
1106619	Small	A7030, A7035
1106620	Medium	A7030, A7035
1106621	Large	A7030, A7035



Available in 4 sizes:



ComfortGel Blue Full



ComfortGel Blue Full full face mask
The design combines from very familiar – and popular – components. Our comfortable forehead cushion and enhanced StabilitySelector combine with our popular blue gel replaceable cushion. The unique sure seal technology (SST) provides a membrane-thin layer over the gel cushion to follow the contours of the face better, allowing wearers to quickly achieve the comfortable seal they need.

Other headgear can be used with this mask, including premium headgear, PS premium headgear, deluxe headgear and small deluxe headgear.

Part number		HCPCS
1040135	ComfortGel Blue Full mask with premium headgear, small	A7030, A7035
1040136	ComfortGel Blue Full mask with premium headgear, medium	A7030, A7035
1040137	ComfortGel Blue Full mask with premium headgear, large	A7030, A7035
1062004*	ComfortGel Blue Full mask with premium headgear, extra-large	A7030, A7035
1002854	ComfortGel Blue Full mask DuoPack with premium headgear, small plus small cushion	A7030, A7035, A7031
1002855	ComfortGel Blue Full mask DuoPack with premium headgear, medium plus medium cushion	A7030, A7035, A7031
1002856	ComfortGel Blue Full mask DuoPack with premium headgear, large plus large cushion	A7030, A7035, A7031
1062005*	ComfortGel Blue Full mask DuoPack with premium headgear, extra-large plus extra large cushion	A7030, A7035, A7031

Accessories		HCPCS
1040568	ComfortGel Blue Full mask without headgear, small	A7030
1040569	ComfortGel Blue Full mask without headgear, medium	A7030
1040570	ComfortGel Blue Full mask without headgear, large	A7030
1062003*	ComfortGel Blue Full mask without headgear, extra-large	A7030
1051032	ComfortGel Blue Full mask without headgear, DuoPack, small plus small cushion	A7030, A7031
1051033	ComfortGel Blue Full mask without headgear, DuoPack, medium plus medium cushion	A7030, A7031
1051034	ComfortGel Blue Full mask without headgear, DuoPack, large plus large cushion	A7030, A7031
1040798	ComfortGel Blue Full cushion, small	A7031
1040799	ComfortGel Blue Full cushion, medium	A7031
1040800	ComfortGel Blue Full cushion, large	A7031
1062030*	ComfortGel Blue Full cushion, extra-large	A7031
1040795	SST flap only, small	A9999
1040796	SST flap only, medium	A9999
1040797	SST flap only, large	A9999
1062031*	SST flap only, extra-large	A9999
1044861	ComfortGel Blue Full cushion and SST flap, small	A7031
1044862	ComfortGel Blue Full cushion and SST flap, medium	A7031
1044863	ComfortGel Blue Full cushion and SST flap, large	A7031
1062032*	ComfortGel Blue Full cushion and SST flap, extra-large	A7031
1040138	ComfortGel Blue Full premium headgear	A7035
1040114	Premium forehead pad	
1049398	Entrainment valve swivel	
1040139	Ball-and-socket Quick Clips, blue, 2 per pk.	
1010872	Port caps, 5 per pk.	
1049507	Sizing gauge	

FitLife

Part number		HCPCS
FitLife with headgear		
1060801	Small	A7030, A7035
1060802	Large	A7030, A7035
1089991	Extra-large	A7030, A7035
FitLife without headgear		
1072530	Small	A7030
1072531	Large	A7030
1089992	Extra-large	A7030
FitLife replacement parts		
1060807	FitLife headgear, small	A7035
1060808	FitLife headgear, large	A7035
1089990	FitLife headgear, extra-large	A7035
1060913	FitLife snap clips, 2 per pk.	
7041	Mask swivel	
1060809	Entrainment elbow	A7045



FitLife
A solution for patients who have difficulty wearing masks or who simply prefer a clear line of sight. FitLife addresses a number of patient issues: Facial hair, dentures, nasal bridge breakdown, facial irregularities and claustrophobia are issues that can greatly impact the ability to establish effective therapy and comfort. FitLife quickly establishes a robust seal around the less pressure-sensitive perimeter of the face, while maintaining an unrestricted view.

Nuance

Nuance gel pillows platform

Philips Respironics, one of the world leaders in gel cushion technology, is bringing the comfort and sealing power of gel to nasal pillows.

As one of the latest in our groundbreaking new line of masks, Nuance gel pillows represent a choice for enhancing sleep therapy for clinicians, home care providers and patients.

Nuance is designed to help patients succeed by providing a comfortable and familiar feel. With the greater comfort and seal of gel pillows, the ability to fit many different size nostrils and choice of a fabric or gel frame, Nuance is truly a great choice in nasal pillows masks.

Key benefits

- Choice between fabric or gel frame for comfort
- Non-slip headgear stays in place throughout the night
- Gel pillows provides excellent comfort and seal
- Lightweight, flexible tubing
- Conforms to different size nostrils



Part number

Nuance		HCPCS
1105160	Nuance gel pillows mask, with headgear	A7034, A7035
1106194	Nuance gel pillows mask, without headgear	A7034
Nuance Pro		
1105167	Nuance Pro gel pillows mask, with headgear	A7034, A7035
1106195	Nuance Pro gel pillows mask, without headgear	A7034
Replacement parts		
1105176	Nuance headgear	A7035
1105405	Nuance RS fabric frame	
1105177	Nuance fabric frame	
1105178	Nuance Pro fabric frame	
1105179	Nuance Pro headgear	A7035
1105179	Nuance Pro frame	
1105173	Nuance/Nuance Pro small gel cushion	A7033
1105174	Nuance/Nuance Pro medium gel cushion	A7033
1105175	Nuance/Nuance Pro large gel cushion	A7033
1105180	Nuance/Nuance Pro swivel tube with exhalation	A7045

102 800.345.6443 www.philips.com/respironics

[Go to table of contents](#)

Nuance. A new choice.

Designed around patients



Nuance fabric frame

- All fabric frame and headgear enhances patient comfort
- Soft fabric experience for better sleep



Nuance Pro Gel frame

- Gel-padded frame and nonslip headgear holds mask in place
- Reduces need for re-adjustment
- High-performance look

Order online, anytime at www.homehealth.estore.philips.com

800.345.6443 103

[Go to table of contents](#)



Accessories/U.S. HCPCS codes

Accessories (cont.)

302175	Chin strap
1012911	Premium chin strap
302425	Deluxe chin strap
1032907	Lightweight CPAP tubing (6 ft.)
7212	Standard CPAP tubing (10 ft.)
1036425	String gauge ring

Codes for nasal masks

A7034	Nasal mask only
A7035	Headgear
A7032	Nasal cushion replacement

Codes for full face masks

A7030	Full face mask only
A7035	Headgear
A7031	Full face cushion replacement

Codes for nasal pillows masks

A7034	Nasal mask only
A7035	Headgear
A7033	Nasal pillow replacement
A7032	Nasal cushion replacement (OptiLife CradleCushion)
A7045	Exhalation port with swivel (GoLife PN 1078407)